

# Calibration and Documentation in Accelerated Weathering and Corrosion Laboratory Testing

Andy Francis – Marketing Director

Bill Tobin – Senior Technical Marketing Specialist

Dave Duecker – Senior Technical Marketing Specialist

Sean Fowler – Senior Technical Director

Q-Lab

[\*View Recorded Presentation\*](#)

# Q-Lab's Standards and Calibration Series

Today is the 1st of a three-part webinar series on standards, calibration, and documentation

All upcoming and archived webinars can be accessed at: [q-lab.com/webinars](https://q-lab.com/webinars)

Date	Topic
01 Sep	Calibration and Documentation
08 Sep	What's New in Standards
15 Sep	How to Run ISO 105-B02

# Administrative Notes

You'll receive a follow-up email from [info@email.q-lab.com](mailto:info@email.q-lab.com) with links to a survey, registration for future webinars, and to download the slides

Use the Q&A feature in Zoom to ask us questions today!



We make testing simple.



**Thank you for attending our webinar!**

We hope you found our webinar on **Calibration and Documentation in Weathering and Corrosion Testing** to be helpful and insightful. The link below will give you access to the slides and recorded webinar.

You can help us continue to provide valuable and high quality content by completing our 3-question [survey](#) about your webinar experience. Every piece of feedback is carefully reviewed by a member of our team.



# Agenda

- Calibration and adjustment
- Benefits of calibration
- Calibration requirements in weathering and corrosion testing
- Q-Lab calibration recommendations
- Calibration documentation
- Accuracy and uncertainty

# Calibration and Adjustment

- **Calibration:** comparison to a known measurement standard. Often traceable back to a national metrological institute (e.g. NIST in the United States).
- **Adjustment:** change made to a measuring system so it outputs proper values of the quantity to be measured.
- *These procedures are often performed together, but not necessarily!*



# Calibration vs Adjustment Example

1. Compare the watch to the official US time, at time.gov. If the watch reads 10:55 and time.gov says 11:00, then your watch has been compared to a known reference and shown to be 5 minutes slow. This is a **calibration**.
2. Fly to a neighboring time zone and move your watch ahead one hour to compensate. This is an **adjustment**. No comparison to a known standard is involved, but this adjustment should get the watch closer to displaying the correct time (unless you got on the wrong flight!)
3. Compare the time on your watch to time.gov and then set your watch to match that time as closely as possible. This is both **calibration and adjustment**.

*Proper weathering testing involves both calibrations and adjustments.*




# Agenda

- Calibration and adjustment
- **Benefits of calibration**
- Calibration requirements in weathering and corrosion testing
- Q-Lab calibration recommendations
- Calibration documentation
- Accuracy and uncertainty

# Calibration Benefits


- Q-Lab and other manufacturers recommend regular calibrations of onboard irradiance and temperature sensors
- Calibration ensures that displayed values of irradiance, temperature, and relative humidity can be trusted
- Even small deviations in measurement and control can lead to significantly different results over a long test period.





# Calibration Benefits

- Calibrating all sensors regularly reduces deviations
- Even with a long test duration, deviation remains small



# The biggest problem with calibrations in weathering testing...

Many end-users neglect to do them, because they often are difficult and/or expensive...


# The biggest improvement we can make for calibrations...

Make it easier and less expensive, so it is more likely that end-users will do them!


# Agenda

- Calibration and adjustment
- Benefits of calibration
- Calibration requirements in weathering and corrosion testing
- Q-Lab calibration recommendations
- Calibration documentation
- Accuracy and uncertainty

# Sensor Terminology



**Onboard (or integrated) sensors** are part of the tester, and are used for measurement and control



**Calibration sensors** are external reference devices used for calibration and adjustment of onboard sensors


# Reference Devices vs. Onboard Sensors

QUV Irradiance

Q-SUN Irradiance

Q-SUN Temp

Onboard sensors



Reference devices



# Calibration in Test Standards

---

- Test standards from ASTM, ISO, and other organizations often require calibrations, but there is not always consistency
- Consider ASTM G155 (xenon arc weathering). This gives very open-ended guidance on calibration:

*6.3 Instrument Calibration*—To ensure standardization and accuracy, the instruments associated with the exposure apparatus (such as timers, thermometers, wet bulb sensors, dry bulb sensors, humidity sensors, UV sensors, radiometers) require periodic calibration to ensure repeatability of test results. Instrument calibration should be traceable to national or international standards. Calibration frequency and procedure should be in accordance with manufacturer's instructions and good laboratory practices.

# Calibration in Test Standards

---

Major, frequently-referenced weathering standards require **at least annual** calibration of onboard sensors. Further guidance is sometimes also provided:

- **ISO 9370** (*Plastics — Instrumental determination of radiant exposure in weathering tests — General guidance and basic test method*) requires onboard sensor calibration checks more frequently than annually, but does not define “frequent.”
  - **ISO 4892-1** (*Plastics — Methods of exposure to laboratory light sources — Part 1: General guidance*) recommends onboard sensor calibration checks more frequently than annually, but does not require them.
  - **ASTM G151** (*Standard Practice for Exposing Materials in Accelerated Test Devices that Use Laboratory Light Sources*) requires onboard sensor calibration checks every time a lamp, optical filter, or test cycle is changed, and annual NMI-traceable calibration and adjustment of reference radiometers.
  - **ISO 9370** and **ASTM G151** both require users to follow manufacturer’s recommendations.
- **It’s not always clear how to meet this guidance!**

# Agenda

- Calibration and adjustment
- Benefits of calibration
- Calibration requirements in weathering and corrosion testing
- **Q-Lab calibration recommendations**
- Calibration documentation
- Accuracy and uncertainty



# Q-Lab's calibration recommendations


- Q-Lab recommends frequent calibrations at regular intervals, for irradiance, temperature, relative humidity, & other parameters
  - We champion this within the standards community
- Technical documentation gives clear guidance on what to do and when
- Q-Lab's **Universal Calibrator** system makes calibration and adjustment simple in QUV and Q-SUN testers

# Universal Calibrator System

## Q-SUN irradiance



## Q-SUN temperature



## QUV irradiance



# Universal Calibrator Benefits

- Practicality
  - Simplest and most cost-effective system on the market
  - Doesn't require an expensive service visit
  - Compatible with **all** Q-Lab testers
- Accuracy
  - Proven, accurate, NIST-traceable calibration method
  - All calibrations of reference devices are performed under test conditions
    - This solidifies the traceability chain
  - All calibrations of onboard sensors performed under actual test condition using reference devices of same type as the onboard sensors

# Calibration: QUV Accelerated Weathering Testers

Sensor	Calibration Frequency	Calibration Instrument
Irradiance	500 hours	UC10 Smart Sensor
Black Panel Temperature	6 months	Calibrated Reference Thermometer

- More frequent calibration is perfectly OK!
- Worth recalibrating also when test cycle, lamps, or sensors are changed



# Calibration: Q-SUN Accelerated Weathering Testers

Sensor	Calibration Frequency	Calibration Instrument
Irradiance	500 hours	UC20 Smart Sensor
Black Panel Temperature	6 months	UC202 Smart Sensor
Chamber Air Temperature Relative Humidity	12 months	<i>Replacement</i>

- More frequent calibration is perfectly OK!
- Worth recalibrating also when test cycle, lamps, filters, or sensors are changed



# Calibration: Q-FOG Corrosion Testers

Sensor	Calibration Frequency	Calibration Instrument
Temperature (Chamber Air, Wet/Dry Bulb)	6 months	Reference thermometer
Shower and Fog deposition	Standards-dependent	Collections devices
Shower flow	6 months	Tester sensors


Q-FOG calibrations, especially for collections, can vary depending on test




# Calibration: Q-FOG Corrosion Verification



Corrosion coupons



Pluviometry



Independent verification

- Widely-used **calibration** techniques (usually with **adjustment**) to ensure that the corrosivity delivered by the tester is what is expected.
- Different style than an onboard sensor calibration, but a calibration nonetheless

# Why Use Q-Lab's Irradiance and Temperature Devices for Calibration?

- Irradiance (UC10, UC20)
  - Calibrated to specific spectrum of xenon lamp/filter combination or UV lamp
  - Eliminates effect of spectral mismatch
  - Avoids saturation issues with off-the-shelf spectrophotometers
- Temperature (UC202)
  - Accounts for radiant heating from xenon lamps
  - Includes convective cooling effects


Reference thermometers and generic UV measurement devices do not calibrate onboard sensors properly for their service environment



# Spectral Mismatch

## Optical Filters in Q-SUN Testers

- Differences between optical filters can lead to mis-calibration if not accounted for
- This is known as spectral mismatch



# Spectral Mismatch

- If a UC20 calibrated for **Daylight-Q** filters measures a tester with **Window-Q** filters installed, an actual irradiance of  $1.00 \text{ W/m}^2/\text{nm}$  @340 nm will measure only  $0.94 \text{ W/m}^2/\text{nm}$
- That spectral mismatch-induced error is 6%
  - This can significantly affect test results!
  - The opposite scenario (i.e. calibrating Daylight-Q with a Window-Q reference device) causes the reverse issue

Imagine errors possible for a device not designed for Q-SUN testers at all!

# Agenda

- Calibration and adjustment
- Benefits of calibration
- Calibration requirements in weathering and corrosion testing
- Q-Lab calibration recommendations
- **Calibration documentation**
- Accuracy and uncertainty

# Q-Lab Global Calibration Locations

## ISO 17025-Accredited


- US
- UK
- China
- Germany



# Calibration Certificates


## Onboard Sensors

Form U-40332-L, Rev. 17, 20 JAN 2021



### Calibration Certificate

**Irradiance Sensors used in QUV and Q-SUN Test Chambers**



The irradiance sensors in this tester are certified to have been calibrated using instrumentation traceable to the National Institute of Standards and Technology (NIST). Irradiance calibration is done in accordance with the following internal Q-Lab procedure: LP-CAL-04. This calibration is ISO 17025 accredited by A2LA.

**Traceability Chain**

QUV or Q-SUN irradiance sensor → radiometer → 3 master radiometers → spectroradiometer → FEL Lamp → NIST

**Customer Name and Address:**

Model:

Serial No:

Prior Calibration Date:

**Working Standard Radiometer:**

Type:

Serial No:

Calibration Date:

**Tester Information:**

Model:

Serial No:

Prior Calibration Date:

**Working Standard Radiometer:**

Type:

Serial No:

Calibration Date:

**Calibration Performed:**

Calibration Date:

Valid Until:

Q-Lab recommends annual, accredited calibrations for irradiance sensors. As such, Q-Lab specifies a default "Valid Until" entry of 12 months from the sensor's original calibration date (ISO 17025, however, allows any customer to request a different or blank date at their complete discretion).

Q-Lab also recommends that the user perform a calibration verification check every 500 light hours. Light hours since the last calibration can be found by pushing the "h" button and setting to 010 on the QUV tester or 05 on Q-SUN testers.

Black Panel temperature:  °C

Chamber air temperature (Q-SUN only):  °C

Lamp type (QUV only):

UV filter type (Q-SUN only):

Irradiance sensor wavelength (Q-SUN only):

	<input type="checkbox"/> Irradiance ( W/m <sup>2</sup> at nm** )				<input type="checkbox"/> Illuminance (lux)			
	Ch. 1	Ch. 2	Ch. 3	Ch. 4	Ch. 1	Ch. 2	Ch. 3	Ch. 4
<b>As Found</b>								
Zero – reference radiometer								
Zero – unit under test*								
Span – reference radiometer								
Span – unit under test								
<b>As Left</b>								
Zero – reference radiometer								
Zero – unit under test*								
Span – reference radiometer								
Span – unit under test								

\* All Q-SUN and QUV tester series 73 and higher or testers with upgraded TEBR. \*\*Irradiance unit for UVC-254 is mW/cm<sup>2</sup>. The calibration uncertainty at a 95% confidence level using a coverage factor of k=2 is ±5.5%.


Calibration Performed By:

Approved By:

This certificate shall not be reproduced except in full without the approval of Q-Lab Corporation. | Page 1 of 1


Q-Lab Corporation | 800 Canterbury Road | Westlake, OH 44145 USA | Repair@q-lab.com | +1-440-835-8700

Form LP-CAL-01, Rev. 20, 20 JAN 2021



### Calibration Certificate

**Temperature and RH Sensors for Use in QUV, QCT, Q-SUN & Q-FOG Test Chambers**



The thermometer (and humidity sensor for Q-SUN Xe-3 Series 1-22) in this tester is certified to have been calibrated using instrumentation traceable to the National Institute of Standards and Technology (NIST). Temperature calibration is done in accordance with internal Q-Lab procedure LP-CAL-01, LP-CAL-01b, and LP-CAL-01c and the temperature scale is ITS-90. Relative humidity calibration (if applicable) is done by direct comparison to a calibrated NIST traceable hygrometer, internal Q-Lab procedure LP-CAL-02. This calibration is ISO 17025 accredited by A2LA.

**Traceability Chains**

temperature sensor in QUV, QCT, Q-SUN (non-BP), Q-FOG and RH sensor in Q-SUN Xe-3 series 1-22 → working standard thermometer (or hygrometer) → NIST

Q-SUN BP temperature sensor → black panel thermometer → master black panel thermometer → reference standard thermometer → NIST

**Customer Name and Address:**

Model:

Serial No:

Prior Calibration Date:

**Tester Information:**

Model:

Serial No:

Prior Calibration Date:

**Calibration Performed:**

Calibration Date:

Valid Until:

Q-Lab recommends annual, accredited calibrations for temperature and RH sensors. As such, Q-Lab specifies a default "Valid Until" entry of 12 months from the sensor's original calibration date. ISO 17025, however, allows any customer to request a different or blank date at their complete discretion.

**Working Standards Used:**

Type of Device	Manufacturer and Model		Serial Number	
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

Sensor Type:	1		2		1		2	
	1	2	1	2	1	2	1	2
<b>Calibration Point:</b>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<b>Reference Thermometer (or Hygrometer Reading):</b>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<b>Tester Temperature (or RH Reading) As Found:</b>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<b>Tester Temperature (or RH Reading) As Left:</b>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

The calibration uncertainties at a 95% confidence level using a coverage factor of k=2 are as follows:  
 Temperature: ± 0.25°C, Q-SUN Black Panel: ± 1.9°C, Relative Humidity: ±3.4%

Calibration Performed By:


Approved By:

This certificate shall not be reproduced except in full without the approval of Q-Lab Corporation. | Page 1 of 1


Q-Lab Corporation | 800 Canterbury Road | Westlake, OH 44145 USA | Repair@q-lab.com | +1-440-835-8700

# Traceability Chain

## Irradiance





Form U-40332-L, Rev. 17, 20 JAN 2021



**Q-LAB**

### Calibration Certificate

Irradiance Sensors used in  
QUV and Q-SUN Test Chambers

The irradiance sensors in this tester are certified to have been calibrated using instrumentation traceable to the National Institute of Standards and Technology (NIST). Irradiance calibration is done in accordance with the following internal Q-Lab procedure: LP-CAL-04. This calibration is ISO 17025, accredited by A2LA.

**Traceability Chain**

QUV or Q-SUN irradiance sensor
⇒
radiometer
⇒
3 master radiometers
⇒
spectroradiometer
⇒
FEL Lamp
⇒
NIST

<p><b>Customer Name and Address:</b></p> <div style="border: 1px solid gray; height: 40px; width: 100%;"></div>	<p><b>Tester Information:</b></p> <p>Model: <input type="text"/></p> <p>Serial No: <input type="text"/></p> <p>Prior Calibration Date: <input type="text"/></p> <p><b>Working Standard Radiometer:</b></p> <p>Type: <input type="text"/></p> <p>Serial No: <input type="text"/></p> <p>Calibration Date: <input type="text"/></p>	<p><b>Calibration Performed:</b> <input type="text"/></p> <p>Calibration Date: <input type="text"/></p> <p>Valid Until: <input type="text"/></p> <p><small>Q-Lab recommends annual, accredited calibrations for irradiance sensors. As such, Q-Lab specifies a default "Valid Until" entry of 12 months from the sensor's original calibration date. ISO 17025, however, allows any customer to request a different or blank date at their complete discretion.</small></p> <p><small>Q-Lab also recommends that the user performs a calibration verification check every 500 light hours. Light hours since the last calibration can be found by pushing the ? button and scrolling to D10 on the QUV tester or D5 on Q-SUN testers.</small></p>
---	---	---

Black Panel temperature:  °C

Lamp type (QUV only):

Chamber air temperature (Q-SUN only):  °C

UV filter type (Q-SUN only):

Irradiance sensor wavelength (Q-SUN only):

As Found	<input type="checkbox"/> Irradiance (W/m <sup>2</sup> at nm)**	Ch. 1	Ch. 2	Ch. 3	Ch. 4
Zero – reference radiometer					
Zero – unit under test*					
Span – reference radiometer					
Span – unit under test					
<b>As Left</b>		Ch. 1	Ch. 2	Ch. 3	Ch. 4
Zero – reference radiometer					
Zero – unit under test*					
Span – reference radiometer					
Span – unit under test					

\* All Q-SUN and QUV tester series 73 and higher or testers with upgraded TEBR. \*\*Irradiance unit for UVC-254 is mW/cm<sup>2</sup>. The calibration uncertainty at a 95% confidence level using a coverage factor of k=2 is ±6.5%.

Calibration Performed By:

Approved By:


This certificate shall not be reproduced except in full without the approval of Q-Lab Corporation. | Page 1 of 1

Q-Lab Corporation | 800 Canterbury Road | Westlake, OH 44145 USA | Repair@q-lab.com | +1-440-835-8700


# Traceability Chain

## Temperature


### Black Panel



### CAT/RH Sensor





Form LF-CAL-01, Rev. 20, 20 JAN 2021



## Calibration Certificate

Temperature and RH Sensors for Use in  
QUV, QCT, Q-SUN & Q-FOG Test Chambers

ACCREDITED  
CERT #038101

The thermometer (and humidity sensor for Q-SUN Xe-3 Series 1-22) in this tester is certified to have been calibrated using instrumentation traceable to the National Institute of Standards and Technology (NIST). Temperature calibration is done in accordance with internal Q-Lab procedure LP-CAL-01, LP-CAL-01b, and LP-CAL-01c and the temperature scale is ITS-90. Relative humidity calibration (if applicable) is done by direct comparison to a calibrated NIST traceable hygrometer, internal Q-Lab procedure LP-CAL-02. This calibration is ISO 17025, accredited by A2LA.

**Traceability Chains**

temperature sensor in QUV, QCT, Q-SUN (non-BP), Q-FOG and RH sensor in Q-SUN Xe-3 series 1-22

→

working standard thermometer (or hygrometer)

→

NIST

Q-SUN BP temperature sensor

→

black panel thermometer

→

master black panel thermometer

→

reference standard thermometer

→

NIST

<b>Customer Name and Address:</b>	<b>Tester Information:</b> Model: Q-FOG GRH600-HTC Serial No: Prior Calibration Date:	<b>Calibration Performed:</b> Calibration Date: Valid Until: <small>Q-Lab recommends annual, accredited calibrations for temperature and RH sensors. As such, Q-Lab specifies a default "Valid Until" entry of 12 months from the sensor's original calibration date. ISO 17025, however, allows any customer to request a different or blank date at their complete discretion.</small>
-----------------------------------	--	---

**Working Standards Used:**

Type of Device	Manufacturer and Model	Serial Number

Sensor Type:	1		2		1		2		1		2	
Calibration Point:												
Reference Thermometer (or Hygrometer Reading):												
Tester Temperature (or RH Reading) As Found:												
Tester Temperature (or RH Reading) As Left:												

The calibration uncertainties at a 95% confidence level using a coverage factor of k=2 are as follows:  
Temperature: ± 0.25°C, Q-SUN Black Panel: ± 1.9°C, Relative Humidity: ±3.4%


<b>Calibration Performed By:</b>	<b>Approved By:</b>

This certificate shall not be reproduced except in full without the approval of Q-Lab Corporation. | Page 1 of 1

Q-Lab Corporation | 800 Canterbury Road | Westlake, OH 44145 USA | Repair@q-lab.com | +1-440-835-8700

# Calibration Temperature Sensors



Issue Date: 30-Aug-2022 form UC-110116-L rev. 5, 21-Aug-2020



**Q-LAB**

Calibration Certificate #030522182992411UC202/BP

UC202 Black Panel Calibration Thermometer for  
Use with Q-SUN Xenon Test Chambers

ACCREDITED  
CERT #2383.01

This Black Panel Calibration Thermometer has been calibrated to test procedure UC-110002-T using instrumentation traceable to the National Institute of Standards and Technology (NIST). It should only be used to calibrate the black panel temperature sensor in Q-SUN Xenon Test Chambers. The Q-Lab Corporation calibration lab is ISO 17025 accredited by A2LA.

Traceability Chain

UC202 Calibration Thermometer → 
 Master Black Panel Thermometer → 
 Reference Standard Thermometer → 
 NIST

Q-Lab Corporation  
Repair Dept  
800 Canterbury Road  
Westlake, OH 44145  
United States

Serial No: 18-29924-11-UC202/BP

Prior Calibration: 20-Jan-2022

Calibration Date: 3-May-2022

Valid Until: 3-May-2023

Certificate will be made valid for 12 months from date of calibration unless otherwise specified by the customer at the time of order placement.

Master Black Panel Thermometer Serial Number: 03-0262-3-BP  
Reference Standard Thermometer Serial Number: 200588533

MASTER DATA	Initial			Final		
	1	2	average	1	2	average
zero (ice bath)	-	-	0.0°C	-	-	0.0°C
full range (in Q-SUN)	89.2°C	92.3°C	90.8°C (A)	88.8°C	91.9°C	90.4°C (C)
mid range (in Q-SUN)	-	-	-	48.5°C	50.4°C	49.5°C

UNIT UNDER TEST DATA	Initial (as received)			Final (after calibration)		
	1	2	average	1	2	average
zero (ice bath)	-	-	0.0°C	-	-	0.0°C
full range (in Q-SUN)	87.6°C	90.3°C	89.0°C (B)	89.0°C	91.6°C	90.3°C (D)
mid range (in Q-SUN)	-	-	-	48.7°C	49.9°C	49.3°C
change ((D-C) - (B-A))	-	-	-	-	-	1.8°C

Condition when received: Good. Proceed with Calibration. The calibration uncertainty at a 95% confidence level using a coverage factor of k=2 is ± 1.5 °C



Comments: R8



# Calibration Certificate: Temperature Sensors

## Traceability Chain

Issue Date: 30-Aug-2022 form UC-110116-L rev. 5, 21-Aug-2020


 Calibration Certificate #030522182992411UC202/BP  
UC202 Black Panel Calibration Thermometer for  
Use with Q-SUN Xenon Test Chambers 

This Black Panel Calibration Thermometer has been calibrated to test procedure UC-110002-T using instrumentation traceable to the National Institute of Standards and Technology (NIST). It should only be used to calibrate the black panel temperature sensor in Q-SUN Xenon Test Chambers. The Q-Lab Corporation calibration lab is ISO 17025 accredited by A2LA.

**UC202 Calibration Thermometer** ⇌ **Master Black Panel Thermometer** ⇌ **Reference Standard Thermometer** ⇌ **NIST**

<b>Q-Lab Corporation</b> Repair Dept 800 Canterbury Road Westlake, OH 44145 United States	Serial No: <u>18-29924-11-UC202/BP</u> Prior Calibration: <u>20-Jan-2022</u>	Calibration Date: <u>3-May-2022</u> Valid Until: <u>3-May-2023</u>  Certificate will be made valid for 12 months from date of calibration unless otherwise specified by the customer at the time of order placement.
---	---	---

Master Black Panel Thermometer Serial Number: 03-0262-3-BP  
Reference Standard Thermometer Serial Number: 200588533



# Calibration Certificate: Temperature Sensors

## Calibration Validity

Q-Lab Corporation  
Repair Dept  
800 Canterbury Road  
Westlake, OH 44145  
United States

Serial No: 18-29924-11-UC202/BP  
Prior Calibration: 20-Jan-2022

Calibration Date: 3-May-2022  
Valid Until: 3-May-2023

Certificate will be made valid for 12 months from date of calibration unless otherwise specified by the customer at the time of order placement.

Master Black Panel Thermometer Serial Number: 03-0262-3-BP  
Reference Standard Thermometer Serial Number: 200588533

# Calibration Certificate: Temperature Sensors

## Test Data

MASTER DATA	Initial			Final		
	1	2	average	1	2	average
zero (ice bath)	-	-	0.0°C	-	-	0.0°C
full range (in Q-SUN)	89.2°C	92.3°C	90.8°C (A)	88.8°C	91.9°C	90.4°C (C)
mid range (in Q-SUN)	-	-	-	48.5°C	50.4°C	49.5°C

UNIT UNDER TEST DATA	Initial (as received)			Final (after calibration)		
	1	2	average	1	2	average
zero (ice bath)	-	-	0.0°C	-	-	0.0°C
full range (in Q-SUN)	87.6°C	90.3°C	89.0°C (B)	89.0°C	91.6°C	90.3°C (D)
mid range (in Q-SUN)	-	-	-	48.7°C	49.9°C	49.3°C
change ((D-C) - (B-A))	-	-	-	-	-	1.8°C

Condition when received: Good. Proceed with Calibration.


Comments:

R8

The calibration uncertainty at a 95% confidence level using a coverage factor of k=2 is  $\pm 1.5$  °C



# Calibration Irradiance Sensors (Radiometers)

Issue Date: 30-Aug-2022 form UC-110100-L rev. 3, 21-Aug-2020



**Q-LAB**

Calibration Certificate #300822202170413UC20/340  
 UC20/340 Radiometer for Use with the SOLAR EYE  
 Irradiance Control System on Q-SUN Xenon Test  
 Chambers

This radiometer is certified to have been calibrated using instrumentation traceable to the National Institute of Standards and Technology (NIST) and in accordance with ASTM test method G130 sec. 6.4. It should only be used to measure xenon lamps in Q-SUN Xenon Test Chambers. Q-Lab's calibration lab is ISO 17025 accredited by A2LA for the calibration of radiometers.

Traceability Chain

```

            graph LR
            A[radiometer] --> B[3 master radiometers]
            B --> C[spectroradiometer]
            C --> D[FEL Lamp]
            D --> E[NIST]
            
```

Q-Lab Corporation Acceptance Testing 800 Canterbury Road Westlake, OH 44145 United States	Serial No: <u>20-21704-13-UC20/340</u> Prior Calibration: <u>26-Aug-2021</u>	Calibration Date: <u>30-Aug-2022</u> Valid Until: <u>30-Aug-2023</u> <small>ASTM G151, sec. 5.1.7.4 and ISO 9370 sec. 6.1.1 require annual calibration. As such, Q-Lab specifies a default "Valid Until" entry of 12 months from the radiometer's original calibration date. ISO 17025, however, allows users to request a different or blank date at their discretion.</small>
---	---	---

Master Radiometer Serial Numbers: #1: 20-19803-12-UC20/340 #2: 20-19809-12-UC20/340 #3: 19-32937-12-UC20/340

MASTER DATA	Irradiance (W/m <sup>2</sup> /nm @340nm)			
	#1	#2	#3	Avg.
Zero	0.00	0.00	0.00	0.00
SPAN, Full-Range	1.30	1.28	1.27	1.28
SPAN, Mid-Range	0.35	0.34	0.34	0.34




UNIT UNDER TEST DATA		
As Received - ZERO		0.00
As Received - SPAN, Full-Range	(A)	1.40
After Cleaning - ZERO		0.00
After Cleaning - SPAN, Full-Range	(B)	1.41
After Calibration - ZERO		0.00
After Calibration - SPAN, Full-Range	(C)	1.28
After Calibration - SPAN, Mid-Range		0.34
% Change After Cleaning - SPAN, Full-Range	(B-A)/A	1%
% Change After Cleaning and Calibration - SPAN, Full-Range (C-A)/A		-9%

Condition when received: Good, Proceed with Calibration. The calibration uncertainty at a  
 Comments: Working Standard R9 95% confidence level using a  
coverage factor of k=2 is ± 5.5%

# Calibration Certificate: Irradiance Sensors

## Traceability Chain


Issue Date: 30-Aug-2022 form UC-110100-L rev. 3, 21-Aug-2020

 Calibration Certificate #300822202170413UC20/340  
UC20/340 Radiometer for Use with the SOLAR EYE  
Irradiance Control System on Q-SUN Xenon Test  
Chambers  

This radiometer is certified to have been calibrated using instrumentation traceable to the National Institute of Standards and Technology (NIST) and in accordance with ASTM test method G130 sec. 6.4. It should only be used to measure xenon lamps in Q-SUN Xenon Test Chambers. Q-Lab's calibration lab is ISO 17025 accredited by A2LA for the calibration of radiometers.

**Traceability Chain**

```
graph LR; A[radiometer] --> B[3 master radiometers]; B --> C[spectroradiometer]; C --> D[FEL Lamp]; D --> E[NIST]
```



# Calibration Certificate: Irradiance Sensors

## Calibration Validity

Q-Lab Corporation  
Acceptance Testing  
800 Canterbury Road  
Westlake, OH 44145  
United States

Serial No: 20-21704-13-UC20/340  
Prior Calibration: 26-Aug-2021

Calibration Date: 30-Aug-2022  
Valid Until: 30-Aug-2023

ASTM G151, sec. 5.1.7.4 and ISO 9370 sec. 6.1.1 require annual calibration. As such, Q-Lab specifies a default "Valid Until" entry of 12 months from the radiometer's original calibration date. ISO 17025, however, allows users to request a different or blank date at their discretion.

Master Radiometer Serial Numbers: #1: 20-19803-12-UC20/340 #2: 20-19809-12-UC20/340 #3: 19-32937-12-UC20/340

# Calibration Certificate: Irradiance Sensors

## Test Data

### Q-SUN UC20

MASTER DATA	Irradiance (W/m <sup>2</sup> /nm @340nm)			
	#1	#2	#3	Avg.
Zero	0.00	0.00	0.00	0.00
SPAN, Full-Range	1.30	1.28	1.27	1.28
SPAN, Mid-Range	0.35	0.34	0.34	0.34
UNIT UNDER TEST DATA				
As Received - ZERO				0.00
As Received - SPAN, Full-Range (A)				1.40
After Cleaning - ZERO				0.00
After Cleaning - SPAN, Full-Range (B)				1.41
After Calibration - ZERO				0.00
After Calibration - SPAN, Full-Range (C)				1.28
After Calibration - SPAN, Mid-Range				0.34
% Change After Cleaning - SPAN, Full-Range (B-A)/A				1%
% Change After Cleaning and Calibration - SPAN, Full-Range (C-A)/A				-9%

### QUV UC10

MASTER DATA (calibration version 2)	Irradiance UV-A Lamps (W/m <sup>2</sup> /nm @340nm)				UV-B Lamps (W/m <sup>2</sup> /nm @310nm)			
	#1	#2	#3	Avg.	#1	#2	#3	Avg.
Zero								
SPAN, Full-Range								
SPAN, Mid-Range								
UNIT UNDER TEST DATA (calibration version 2)								
As Received - ZERO								
As Received - SPAN, Full-Range (A)								
After Cleaning - ZERO								
After Cleaning - SPAN, Full-Range (B)								
After Calibration - ZERO								
After Calibration - SPAN, Full-Range (C)								
After Calibration - SPAN, Mid-Range								

Condition when received:  
Comments:

Good. Proceed with Calibration.  
Working Standard R9

The calibration uncertainty at a  
95% confidence level using a  
coverage factor of k=2 is ± 5.5%

# Agenda

- Calibration and adjustment
- Benefits of calibration
- Calibration requirements in weathering and corrosion testing
- Q-Lab calibration recommendations
- Calibration documentation
- Accuracy and uncertainty



# Calibration Uncertainty Budget

- The accuracy of onboard sensors is determined by the accuracy of the references used to calibrate them.
- This is calculated by preparing an uncertainty budget, shown here for UC20 irradiance smart sensors

UC 20 source of uncertainty	Est. uncertainty (%) 2 Std-Dev	Std. uncertainty (%) 1 Std-Dev
Parameter #1	1.16	0.58
Parameter #2	1.00	0.50
Parameter #3	1.00	0.50
Parameter #4	0.08	0.04
Parameter #5	0.08	0.04
Parameter #6	0.20	0.10
Parameter #7	2.50	1.25
Parameter #8	1.00	0.58
Parameter #9	0.28	0.14
Parameter #10	4.80	2.40
Parameter #11	0.30	0.20
Parameter #12	0.60	0.30
Parameter #13	1.43	0.71
Parameter #14	2.00	1.00
Parameter #15	0.70	0.35
Parameter #16	2.86	1.43
Parameter #17	1.43	0.71
<b>Combined Uncertainty</b>	<b>±5.6%</b>	<b>±2.8%</b>

# Accuracy and Uncertainty

- A very common question in calibration: *What is the accuracy of my sensors?*
- Consider BP/IBP sensors:
  - Calibrated with a  $k=2$  measurement uncertainty of  $\pm 1.5$  °C. This is stated on the device calibration certificate.
  - Calibrated under a radiant load from a xenon-arc lamp and are intended for use specifically in a Q-SUN Xenon Arc tester. This is the main source of uncertainty.
  - We calibrate these devices using reference black panel sensors, and our laboratory accredited to perform the calibration in accordance with ISO 17025.

The accuracy of onboard sensors is determined by the accuracy of the references used to calibrate them

# Summary – Calibration and Documentation

- Calibration of onboard sensors is critical for correct, consistent performance of weathering and corrosion testers
- Calibration requirements from test standards are not always very descriptive
- The Universal Calibrator system makes irradiance and temperature calibration simple for QUV and Q-SUN testers
- Extensive documentation is provided with calibration certificates for onboard sensors and calibration devices, including:
  - Traceability chain
  - Customer information and calibration validity
  - Test conditions
  - Reference devices and test data
- Sensor accuracy is determined by the accuracy of the devices used to calibrate them

# Thank you for your attention!

## Questions?

Send your inquiry to:  
[info@q-lab.com](mailto:info@q-lab.com)