

Q-SUN Xenon Arc Testers

Operation and Maintenance

Andy Francis – Marketing Director

Bill Tobin – Senior Technical Marketing Specialist

Dave Duecker – Senior Technical Marketing Specialist

Sean Fowler – Senior Technical Director

Q-Lab



[View Recorded Presentation](#)

Q-Lab's Operator Training Series

Today is the first of a three-part webinar series on basic operation of our weathering and corrosion testers


All upcoming and archived webinars can be accessed at:
q-lab.com/webinars

Date	Topic
06 Oct	QUV
13 Oct	Q-SUN
20 Oct	Q-FOG


Administrative Notes

You'll receive a follow-up email from info@email.q-lab.com with links to a survey, registration for future webinars, and to download the slides

Use the Q&A feature in Zoom to ask us questions today!



We make testing simple.



Thank you for attending our webinar!

We hope you found our webinar on *Q-SUN Xenon Arc Tester Operation and Maintenance* to be helpful and insightful. The link below will give you access to the slides and recorded webinar.

You can help us continue to provide valuable and high quality content by completing our [3-question survey](#) about your webinar experience. Every piece of feedback is carefully reviewed by a member of our team. In fact, today's webinar was created as a direct result of customer feedback from previous webinar surveys!

We consistently hold seminars and webinars about weathering, corrosion, standards and more. The best way to keep up with news and events is by following us on [Facebook](#), [Twitter](#) and [LinkedIn](#).

Q-SUN Xenon Arc Models


Simulate light (outdoor direct sunlight, indoor filtered light); heat (elevated temperature), and water (humidity and water spray)



Q-SUN Topics

- Safety
- Features and Functions
- Running a Test
- Calibration
- Maintenance

Xenon Lamps



Xenon Lamps give off UV light and should not be viewed directly

Interlocks stop tester when door opened

Interlock Switches



Xe-1



Xe-2
(overhead view)



Xe-3

Q-SUN Hazards



Electrical Shock



Hot Surface



Topics

- Safety
- Features and Functions
- Running a Test
- Calibration
- Maintenance

Q-SUN Overview (Flat Array)

- 1) User interface
- 2) USB port for data transfer
- 3) Xenon lamps with irradiance control
- 4) Optical filters
- 5) Water spray
- 6) Onboard irradiance sensors
- 7) Black Panel Temp sensor
- 8) Specimen holders
- 9) Relative Humidity/CAT sensor



Q-SUN Overview (Rotating Rack)

- 1) User interface
- 2) USB port for data transfer
- 3) Xenon lamps with irradiance control
- 4) Optical filters
- 5) Water spray
- 6) Onboard irradiance sensors
- 7) Black Panel Temp sensor
- 8) Specimen holders
- 9) Relative Humidity/CAT sensor



Q-SUN Optional Features


Q-SUN Model/ Configuration		Q-SUN Tester Features					
		Gen 4 + High Irrad (-E) ⁴	Humidity Control (-H)	Chamber Air Chiller (-C)	Water Spray (-S)	Back Spray (-BS)	Dual Spray (-DS)
Xe-1	Xe-1-B ³						
	Xe-1-BCE	●		●			
	Xe-1-SE	●			●		
	Xe-1-SCE	●		●	●		
	Xe-1-WE	●			● ⁵		
Xe-2	Xe-2-HE	●	●				
	Xe-2-HSE	●	●		●		
	Xe-2-HBSE	●	●		●	●	
Xe-3	Xe-3-H ³		●				
	Xe-3-HCE	●	●	●			
	Xe-3-HSE	●	●		●		
	Xe-3-HSCE	●	●	●	●		
	Xe-3-HBSE	●	●		●	●	
	Xe-3-HBSCE	●	●	●	●	●	
	Xe-3-HDSE	●	●		●		●
	Xe-3-HDSCE	●	●	●	●		●
	Xe-3-HDSBSE	●	●		●	●	●
	Xe-3-HDSBSCE	●	●	●	●	●	●

Q-SUN Light Delivery

Xenon arc lamps



Optical filters





Xe-2 lantern



Xe-1 / Xe-3 flat filter

Optical Filters and Lamps



Xe-1
and
Xe-3




Xe-2



Filter
Lantern

Xenon Arc Spectrum



SOLAR EYE Irradiance Control System


- SOLAR EYE Irradiance Control maintains the same light output at all times
- Ballasts control lamp output
- Allows for excellent repeatability and reproducibility




Onboard Irradiance Sensors



Xe-1



Xe-2
(overhead)




Xe-3

Xe-1 Onboard Temperature Sensors

- Xe-1 Black Panel (BP) or Insulated Black Panel (IBP) Thermometer
- Optional Chamber Air Thermometer (CAT) in Xe-1
- Tester can control temperature by BP or CAT; the other is simply monitored

Black Panel



Chamber Air




Xe-2 Onboard Temperature & RH Sensors

- Xe-2 Black Panel (BP) or Insulated Black Panel (IBP) Thermometer
- One sensor monitors Chamber Air Temperature and Relative Humidity (CAT/RH Sensor), standard in Xe-2
- Tester simultaneously controls BP, CAT, & RH.

CAT / RH
Sensor

IBP
Sensor



Xe-3 Onboard Temperature & RH Sensors

- Xe-3 Black Panel (BP) or Insulated Black Panel (IBP) Thermometer
- One sensor monitors Chamber Air Temperature and Relative Humidity (CAT/RH Sensor), standard in Xe-3
- Tester simultaneously controls BP, CAT, & RH.


BP
Sensor



CAT / RH
Sensor



Temperature Control



Optional Chiller (Xe-1 / Xe-3)



Reduces minimum BPT by:

~10 °C (Xe-3)

~20 °C (Xe-1)

Q-SUN Water Delivery


- Water Spray (optional)
 - Front
 - Back (Xe-2 / Xe-3)
 - Dual (Auxiliary) (Xe-3)
 - Immersion (Xe-1)
- Relative Humidity control (Xe-2 / Xe-3)

Q-SUN Water Quality

- Purified (RO/DI) water is **required**
- Spray water requires even higher purity and also low silica
- Q-Lab follows ASTM G151 recommendations

Water Delivery	Resistivity ($\Omega \cdot \text{cm}$)	Conductivity ($\mu\text{S}/\text{cm}$)	Silica (ppm)	Total Dissolved Solids (ppm)	pH
Spray	> 5 M	< 0.2	< 0.1	< 0.1	6-8
Humidity	> 200 k	< 5.0	No requirement	< 2.5	6-8

Spray Systems (Optional)



Front (Top)




Back (Bottom)

Xe-1 / Xe-3 Water Spray System


- Features
 - Pulse Rate Control
 - Automatic Fault Detector
- Two nozzles, used for:
 - Mist
 - Thermal Shock
 - Erosion



Xe-1



Xe-3



Xe-2 Water Spray System

- Features
 - Pulse Rate Control
 - Automatic Fault Detector
- One nozzle, used for:
 - Mist
 - Thermal Shock
- Second Nozzle for optional Back Spray



Specialized



Xe-1 Immersion



Xe-3 Dual Spray

Relative Humidity Control (Xe-2 / Xe-3)

- Feedback Loop System
- RH/CAT Sensor
- Main Controller
- Humidity Generator
 - Nebulizer in Xe-2
 - Boiler in Xe-3




Xe-2



Xe-3




Q-SUN Specimen Exposure Areas



Xe-1



Xe-2



Xe-3


Q-SUN Topics

- Safety
- Functions of the Tester
- Running a Test
- Calibration
- Maintenance

Step 1: Select an Optical Filter

- Check your test method
 - Almost all test methods describe the filter type
- If you don't test to a specific method:
 - Outdoor applications: Daylight - Q
 - Indoor applications: Window - Q

Q-SUN Optical Filters



Step 2: Select Irradiance Settings



- 340 nm Control Point
 - For outdoor service environments
- 420 nm Control Point
 - For indoor service environments
- TUV Control Point
 - 300-400 nm (wideband, general service environments)
 - Typically used for European standards

Select Irradiance Values

	Xe-1 & Xe-3 Irradiance Values Typical (& Maximum) ^{A,B,C}			Xe-2 Irradiance Values Typical (& Maximum) ^{A,B,C}		
	W/m ² /nm @ 340 nm	W/m ² /nm @ 420 nm	W/m ² @ TUV (300-400 nm)	W/m ² /nm @ 340 nm	W/m ² /nm @ 420 nm	W/m ² @ TUV (300-400 nm)
Daylight-F	0.80 (1.30)	1.50 (2.40)	75 (125)	0.80 (0.95)	1.50 (1.70)	75 (85)
Daylight Q	0.68 (1.10)			0.68 (0.80)		
Extended UV (-Q/B, -Quartz ^D)				0.51 (0.61) ^E		
Daylight-B/B	0.55 (0.85)		70 (108)	0.55 (0.65)		70 (80)
Window (-Q, -B/SL)	-		42 (68)	-		42 (62)
Window (-SF5, -IR, -B04 ^F)						

Achievable irradiance values vary by optical filter, tester type, and control point

Step 3: Select Black Panel Type

Panel	Construction	ASTM Designation	ISO Designation	Temp Range (°C)
 A photograph showing a black rectangular panel with a black handle and a blue pen resting on it. A small metal component is attached to the top left corner.	Black painted stainless steel	Uninsulated Black Panel	Black Panel	45-110
 A photograph showing a black rectangular panel mounted on a white base. A blue pen and a small metal component are resting on the panel.	Black painted stainless steel mounted on 0.6 cm white PVDF	Insulated Black Panel	Black Standard	50-120

Black Panel Temperature Control

- BP/IBP is typically mounted in holder (sometimes directly on tray)
- Irradiance, color, and thickness all affect specimen temperature



Step 4: Programming

- Light
- Dark
- Light + Spray
- Dark + Spray
- Dark + Spray Front and Back
- Dual Spray
- Light + Dual Spray
- Light + Immersion
- Dark + Immersion

Q-SUN Front Control Panel




Dual, Full-Color Touchscreen Displays

USB

LED Status Indicator


Status and Menu Screens



Gen 4 Q-SUN Xe-3
Step 1 Light


	Irradiance (W/m ² @ 340 nm)			Temperature (°C)		RH
	Lamp 1	Lamp 2	Lamp 3	BP/BP	Chamber	(%)
Actual	1.10	1.10	1.10	70	47	59
Set	1.10	1.10	1.10	70	47	--

	Step Time	Test Time	Test Energy	Total Time
	(Hrs:Mins)	(Hrs:Mins)	(KJ:/m ²)	(Hrs)
Elapsed	13:25	61:27	154.1	62
Set	24:00	1000:00	--	








Status screen

Setpoint and actual controls
Test timers



Main Menu


-  Manage Test Timers
-  Settings
-  Manage Cycles
-  Diagnostics
-  Calibrate
-  Contact Q-Lab




Menu screen

Cycle management
Calibration
Settings

Programming Tests



Test Duration



Managing Cycles

Step 5: Calibration

More on this later!


Step 6: Specimen Mounting

- Specimen Capacity
 - Xe-1: 17 (51 × 102 mm)
 - Xe-2: 31 (45 × 132 mm)
 - Xe-3: 55 (51 × 102 mm)
- Type of holder
 - Open-Backed (thick, rigid specimens)
 - Solid-Backed (flexible specimens)
 - Three-dimensional
- Specimen Tray
 - Solid
 - Mesh (open)
- Masking
 - Commonly used in textile testing

Specimen Holders Xe-1 / Xe-3




Panel Holders



3D Specimen Holders

Specimen Mounting Trays Xe-3




Solid Tray




Open Mesh Tray


Q-SUN Xe-2 Specimen Holders



Open-backed




Close-backed



Textile masks

Mounting Flexible Specimens



Textiles



Thin Films

Step 7: Running the Test







- Specimen Repositioning
- Tester status and monitoring
- Data Logging via VIRTUAL STRIPCHART (optional)
- Calibration

Specimen Repositioning

- Ensures best repeatability and reproducibility
- Perform at least 4 times per test (weekly for long tests)
- Important for both rotating rack and flat array testers



LED Status Indicator

Color	Appearance	Meaning
Red		<i>Flashing</i> Error, test stopped
Yellow		<i>Flashing</i> Notification, test still running
White		Static Power on, stopped, no active error
Green		Static Test running, no active error
Blue		<i>Flashing</i> Test completed
Magenta		<i>Flashing</i> Software install or VSC transfer

Tester Status

Running Cycle A: ASTM G154 C...
Step 1 UV

	Irradiance (W/m ² @ 340 nm)	Temperature (°C)		RH (%)
	Lamp 1	BP/IBP	Chamber	
Actual	1.08	50	39	40
Set	1.10	50	39	40

	Step Time (Hrs.Mins)	Test Time (Hrs.Mins)	Test Energy (KJ/m ²)	Total Time (Hrs)
	Elapsed	9:26	154:01	609.3
Set	24:00	1000:00	-	-

STOP

Test Completed A: ASTM G154 C...
Step 1 UV

	Irradiance (W/m ² @ 340 nm)	Temperature (°C)	
	Lamp 1	BP/IBP	Chamber
Actual	0.00	25	25
Set	0.35	25	25

	Step Time (Hrs.Mins)	Test Time (Hrs.Mins)	Test Energy (KJ/m ²)	Total Time (Hrs)
	Elapsed	8:08	130:30	164.2
Set	24:00	1000:00	-	-

RUN

Test Stopped A: ASTM G154 C...
Step 1 UV

	Irradiance (W/m ² @ 340 nm)	Temperature (°C)	
	Lamp 1	BP/IBP	Chamber
Actual	0.00	25	25
Set	0.35	25	25

	Step Time (Hrs.Mins)	Test Time (Hrs.Mins)	Test Energy (KJ/m ²)	Total Time (Hrs)
	Elapsed	8:08	130:30	164.2
Set	24:00	1000:00	-	-

RUN



Test Stopped A: Gen 4 Q-SUN Xe-2
Step 1 Light

	Irradiance (W/m ² @ 340 nm)	Temperature (°C)		RH (%)
	Lamp 1	BP/IBP	Chamber	
Actual	0.10	50	39	40
Set	1.10	50	39	40

	Step Time (Hrs.Mins)	Test Time (Hrs.Mins)	Test Energy (KJ/m ²)	Total Time (Hrs)
	Elapsed	9:26	154:01	609.4
Set	24:00	1000:00	-	-

RUN

Notifications and Diagnostics



Q-SUN Topics

- Safety
- Functions of the Tester
- Running a Test
- Calibration
- Maintenance

Q-SUN Calibration

- Irradiance sensor
 - Every 500 light hours
- Black Panel temp sensor
 - Every 6 months
- Chamber air temp sensor
 - Every 12 months
- Calibration best practices
 - Calibrate whenever test conditions are changed
 - Calibrate at operating conditions
 - Always calibrate the lamps before calibrating the Black Panel

Q-SUN Smart Sensor Options



UC20/340



UC202/BP




UC202/IBP

Description	Name	Color Code	Used For
Irradiance Smart Sensor	UC20/340	Green	Calibration Radiometer for 340nm on board sensor
	UC20/420	Red	Calibration Radiometer for 420nm on board sensor
	UC20/TUV	Grey	Calibration Radiometer for TUV on board sensor
	UC20/LUX	Olive	Calibration Radiometer for LUX on board sensor
Temperature Smart Sensor	UC202/BP	Black	Calibration of onboard Black Panel sensor
	UC202/IBP		Calibration of onboard Insulated Black Panel sensor

Irradiance Calibration Ports



Xe-1



Xe-2





Xe-3

Irradiance Calibration



Plug Smart Sensor directly into tester USB port



Select optical filter and *Calibrate*


Black Panel Temp Sensor Calibration

Tester black panels and calibrated black panel are placed adjacent to each other in the chamber



Onboard Sensor in Panel Holder


Onboard Black Panel
Calibration Black Panel



Onboard Sensor on Tray

Black Panel Thermometer Calibration

Let temperature stabilize and press
Calibrate




RH/CAT Sensor Annual Replacement

Xe-2



Xe-3



Xe-1 Chamber Air Temp Sensor Calibration

Using an independent reference device, calibrate the chamber air temperature in hot water bath



Yearly Reference Device Calibration

Irradiance & Black Panel Calibration Devices

- UC20 Irradiance Smart Sensors (replace/recalibrate)
- UC202 Temperature Smart Sensors (replace/recalibrate)

Chamber Air Temp / RH sensor (replace)

Q-SUN Topics

- Safety
- Functions of the Tester
- Running a Test
- Calibration
- Maintenance

Xenon Lamp Replacement Schedule

- Xenon lamps need to be replaced because they “age” (spectral shift to less UV)
- Q-SUN Optical Filters do not age

Q-SUN tester	Irradiance	Warranted Lamp Life (hrs)
Legacy	Typical	1500
Q-SUN “E” Model	Typical	3000
Q-SUN “E” Model	Maximum	1000

Q-SUN Xe-1 / Xe-3 Lamp Replacement

(1) Remove lamp housing from tester



(2) Remove lamp and clean housing




(3) Clean UV filters



Q-SUN Xe-1 / Xe-3 Lamp Replacement


(4)
Replace
lamp



(5) Verify
trigger
finger
contact



(6) Final
lamp
check



(7) Install
replacement
lamp




Q-SUN Xe-2 Lamp Replacement

(1) Remove old lamp



(2) Remove and clean filter lantern with ammonia





Q-SUN Xe-2 Lamp Replacement

(3) Install replacement lamp



(4) Verify trigger finger contact



Routine Maintenance


- Clean / Replace Air Filters (Monthly)
- Clean Humidifier
- Inspect Water Filter
- Clean Spray Nozzles
- Inspect Chamber Wall Reflectors

Every
6 months

Clean or Replace Air Filters



Disposable



Washable (Preferred)

Clean or Replace Air Filters




Xe-1 Air Filter



Chiller Air Filter




Xe-3
Chamber
Air Filter





Xe-3
Ballast
Blower
Air Filter

Clean or Replace Air Filters




Xe-2
Chamber
Air Filter


Xe-2 Lamp
Cooling
Air Filter



Clean Xe-2 Humidifier



Water pan




Reservoir

- Disassemble humidifier water pan and reservoir
- Clean with alcohol or mild detergent

Inspect Water Filter



Xe-3



Xe-2

Clean Spray Nozzles

- Remove and disassemble for cleaning
- Clean in ultra-sonic cleaner or ...
- Clean with anti-scaling cleaner (to remove calcium and magnesium deposits) or...
- Thorough washing and rinsing in detergent



Inspect Chamber Reflectors



Clean - OK



Cloudy - Replace

Clean slip ring track (Xe-2)



- Wipe down with scouring pads
- Clean with isopropyl alcohol



Q-SUN Xe-2 now has a
contactless system for
BP Temp Sensor

Field Calibration Audits, Tester Commissioning, and Customer Education

- This presentation was a condensed version of our Q-SUN operator training. Typical training includes hands-on sessions and further in-depth review of tester components not covered here
- Q-Lab Repair team offers tester audits and field calibrations, in addition to their on-site repair visits and troubleshooting services.
- Q-Lab offers customized training plans that can be catered to each customer's demand and may include *Tester Commissioning* as well as *Weathering 101* and *Atmospheric Corrosion* education
- Contact info@q-lab.com or Repair@q-lab.com for more info

Thank you for your attention!

Questions?

Send your inquiry to:
info@q-lab.com